

Otoport OAE+ABR

for infant screening



The Perfect OAE+ABR Combination

- Otodynamics ABR screening technology
- Gold standard OAE screening and analysis
- Upgrade existing Otoport with the ABR sleeve*

* Ask your dealer for details

 **Otodynamics**
Audiology Systems



Features

The Otodynamics Otoport OAE+ABR Lite is a new handheld screening instrument offering a combination of in-depth ABR and either DPOAE or TEOAE screening. The instrument comprises an Otoport Lite in an ABR sleeve. The ABR sleeve can be purchased separately and fitted retrospectively to upgrade existing Otoports.*

ABR Functions

- Short Chirp stimulation via OAE at 55/s (approx.)
- Screening Stimulus levels 30 to 45dBnHL
- In-the-ear stimulus calibration
- Pre-test physiological and powerline noise analysis
- Helpful on-screen guidance based on test conditions
- Individual three electrode impedance analysis
- Weighted averaging for noise management
- Intelligent powerline interference rejection
- Automated ABR identification using Fsp
- ABR confirmation with configurable template correlation
- ABR waveform and test history for audit trail

OAE Functions

- Automated OAE screening
- Automatic probe fit and noise analysis
- Configurable pass criteria
- SNR on 2, 3 or 4 frequencies
- DPOAE model 2-6kHz 4 frequencies
- TEOAE model 1-4kHz 5 half-octave bands

General Functions

- Light, portable with long battery life
- Easy to use intuitive graphic display and controls
- Customisable testing protocols
- 500 test database storing name, date of birth, gender and notes
- Work list facility
- Otolink PC software with response and waveform viewer
- Data backup to PC
- Facility for user login with password protection
- Admin and User access levels
- Battery Charge indicator and condition gauge
- Recharge via charger or PC
- Stimulus and noise indicator lights
- Hot keys, cursor keys and full alphanumeric keyboard
- Configurable display, backlight, contrast and sound
- Simple or detailed print options

Quality Checks

- Probe QA checks
- Self test of ABR and OAE electronics
- QA test history viewing
- Electrode cable continuity test

Options

- Wireless connectivity
- Linear barcode
- RFID scanner

*ABR upgrade for Colour Otoport available now.
Check with your dealer for availability for black & white model

Otodynamics reserves the right to change specifications and appearances as part of its commitment to product improvement



Screening with Confidence

1 SELECT TEST TYPE

Either ABR or OAE automated screening can be selected through the intuitive menu options.

2 INFORMATIVE SCREENS

The pre-Test screen monitors powerline interference and infant noise. Individual electrodes are highlighted if their impedance is too high, or if their balance is poor

3 'IN THE EAR' CALIBRATION

Every test starts with automated 'in the ear' calibration ensuring accurate results

4 CLEAR TEST PROGRESSION

Rich data and open decision process provides confidence in ABR measurements

5 FAST, AUTOMATED RESULTS

Unique combination of Fsp and Template analysis securely identifies genuine ABRs providing clear, reliable and trusted results



Otolink

Data Management Suite

Downloading result wirelessly or via USB the Otolink suite manages your data easily and efficiently

Test details can be quickly reviewed, printed and archived.

Data can be exported. For example to Oz, Hi Track or Noah in the Americas.

Printing directly to PDF files is available

The Otoport suite also provides patient work list upload facility



